

# Laser Hair Removal Device Technical Comparison

## FDA 510(k) Device Analysis: GEX to OHT Transition

Table 1: Comprehensive Device Comparison with Energy Analysis

Device	510(k) Number	Year	Product Code	User Type	Wavelength (nm)	Beam Size	Area (cm <sup>2</sup> )	Energy Density Range (J/cm <sup>2</sup> )	Max Energy Density (J/cm <sup>2</sup> )	Max Total Energy (J)*
TRIA	K090820	2009	GEX	Consumer	~800	Φ 10mm	0.785	6, 10, 14, 18, 22	22	17.3
SILKPRO	K142845	2015	OHT	Consumer	810	9×9mm	0.81	5, 10, 15, 20, 25	25	20.3
IVYLASER	K230025	2023	OHT	Consumer	808	9×18mm	1.62	1.85, 3.09, 6.17, 9.26, 15.43	15.43	25.0
Diode Laser (Lotuxs V8)	K230090	2023	OHT	Consumer	810	30×10mm	3.0	2.2, 3.4, 4.5, 5.6, 7.0, 9.0	9.0	27.0
Diode Laser (Lotuxs V4)	K230090	2023	OHT	Consumer	810	30×10mm	3.0	2.0, 2.8, 3.5, 4.3, 5.3, 7.0	7.0	21.0

\*Max Total Energy = Max Energy Density × Beam Area

### Key Technical Observations:

**Energy Density Evolution:** The transition from GEX to OHT did not result in universally lower energy densities. The 2015 SILKPRO (OHT) actually had higher maximum energy density (25 J/cm<sup>2</sup>) than the 2009 TRIA (GEX, 22 J/cm<sup>2</sup>). However, later OHT devices show more conservative energy limits.

**Professional vs Consumer Applications:** The 2018 professional GEX device (K181019) maintains the highest energy density range (5-40 J/cm<sup>2</sup>), reflecting its clinical use environment with trained operators.

**Total Energy Output Trends:** Despite similar energy densities, total energy output varies significantly based on beam size. Larger beam sizes in recent OHT devices (like the 3.0 cm<sup>2</sup> Lotuxs devices) can deliver substantial total energy even at moderate energy densities.

**Regulatory Timeline:** The precise transition point was July 1, 2015, with the SILKPRO device serving as the regulatory bridge from GEX to OHT classification for consumer laser hair removal devices.

Device	510(k) Number	Year	Product Code	User Type	Wavelength (nm)	Beam Size	Area (cm <sup>2</sup> )	Energy Density Range (J/cm <sup>2</sup> )	Max Energy Density (J/cm <sup>2</sup> )	Max Total Energy (J)*
Diode Laser System	K181019	2018	GEX	Professional	808	-	1.20	5-40	40	48.0

**Data Sources:** FDA 510(k) Premarket Notification Database **Product Codes:**

- GEX: Powered Laser Surgical Instrument
- OHT: Light Based Over-The-Counter Hair Removal Device